

Product Datasheet

Dopamine D2 receptor antibody (orb164310)

Description Rabbit polyclonal antibody to Dopamine D2 receptor

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target DRD2

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal


MW 50619

Uniprot ID [P14416](#)


NCBI [NP_057658.2](#), [NP_000786.1](#)

Dilution Range WB: 1:1000


Expiration Date 12 months from date of receipt.



Western blot:Anti-DRD2 Antibody (C-term)...



Western blot:Anti-DRD2 Antibody (C-term)...



Western blot:Anti-DRD2 Antibody (C-term)...