

Product Datasheet

ASB3 antibody (orb164319)

Description	Rabbit polyclonal antibody to ASB3 (N-term).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ASB3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57745
Uniprot ID	Q9Y575
NCBI	NP_001157637.1 , NP_057199.1 , NP_665862.1 , NP_001188894.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.