

Product Datasheet

PLRG1 antibody (orb164351)

Description	Rabbit polyclonal antibody to PLRG1 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PLRG1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57194
Uniprot ID	O43660
NCBI	NP_002660.1 , NP_001188493.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.