

Product Datasheet

SRPX antibody (orb164372)

Description	Rabbit polyclonal antibody to SRPX
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SRPX
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51572
Uniprot ID	P78539
NCBI	NP_001164222.1 , NP_006298.1 , NP_001164221.1 , NP_001164223.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.