

Product Datasheet

TTC8 antibody (orb164381)

Description	Rabbit polyclonal antibody to TTC8 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TTC8
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	61534
Uniprot ID	Q8TAM2
NCBI	NP_653197.2 , NP_938052.1 , NP_938051.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.