

Product Datasheet

TTC34 antibody (orb164385)

Description	Rabbit polyclonal antibody to TTC34 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TTC34
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60903
Uniprot ID	A8MYJ7
NCBI	NP_001229601.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.