

Product Datasheet

39692 antibody (orb164389)

Description	Rabbit polyclonal antibody to 39692 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SEPTIN8
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55756
Uniprot ID	Q92599
NCBI	NP_001092281.1 , NP_001092282.1 , NP_001092283.1 , NP_055961.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.