

Product Datasheet

CLCN3 antibody (orb164420)

Description	Rabbit polyclonal antibody to CLCN3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CLCN3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	90966
Uniprot ID	P51790
NCBI	NP_776297.2, NP_001230303.1, NP_001820.2, NP_001230301.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.