

Product Datasheet

CKS1B antibody (orb164422)

Description	Rabbit polyclonal antibody to CKS1B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CKS1B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	9660
Uniprot ID	P61024
NCBI	NP_001817.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.