

Product Datasheet

GNMB antibody (orb164487)

Description	Rabbit polyclonal antibody to GNMB (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GNMB
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	63923
Uniprot ID	Q14956
NCBI	NP_002501.1 , NP_001005340.1
Dilution Range	WB: 1:500
Expiration Date	12 months from date of receipt.