

Product Datasheet

EPDR1 antibody (orb164512)

Description	Rabbit polyclonal antibody to EPDR1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	EPDR1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	25437
Uniprot ID	Q9UM22
NCBI	NP_001229875.2 , NP_001229877.1 , NP_060019.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.