

Product Datasheet

AP4B1 antibody (orb164527)

Description	Rabbit polyclonal antibody to AP4B1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	AP4B1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	83260
Uniprot ID	Q9Y6B7
NCBI	NP_006585.2 , NP_001240782.1 , NP_001240781.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.