

Product Datasheet

PLEKHG5 antibody (orb164532)

Description	Rabbit polyclonal antibody to PLEKHG5 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PLEKHG5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	111231
Uniprot ID	O94827
NCBI	NP_001252521.1 , NP_941374.2 , NP_065682.2 , NP_001252522.1 , NP_001036129.1 , NP_001036130.1 , NP_001252523.1 , NP_001036128.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.