

Product Datasheet

TGF beta Receptor 2 antibody (orb164564)

Description	Rabbit polyclonal antibody to TGF beta Receptor 2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Tgfr2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	64219
Uniprot ID	Q62312
NCBI	NP_033397.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.