

Product Datasheet

PRDM10 antibody (orb164629)

Description	Rabbit polyclonal antibody to PRDM10
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PRDM10
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	130136
Uniprot ID	Q9NQV6
NCBI	NP_955469.1 , NP_955470.1 , NP_064613.2 , NP_955471.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.