

Product Datasheet

Recombinant Human Type-1 angiotensin II receptor-associated protein(AGTRAP) (orb1646346)

Description	Recombinant Human Type-1 angiotensin II receptor-associated protein(AGTRAP)
Species/Host	Bacteria
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.