

Product Datasheet

HIST1H2AB antibody (orb164670)

Description	Rabbit polyclonal antibody to HIST1H2AB (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	H2AC4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	14135
Uniprot ID	P04908
NCBI	NP_066390.1 , NP_003504.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.