

Product Datasheet

MTMR14 antibody (orb164691)

Description	Rabbit polyclonal antibody to MTMR14 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MTMR14
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	72203
Uniprot ID	Q8NCE2
NCBI	NP_001070993.1 , NP_001070994.1 , NP_071930.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.