

Product Datasheet

ZNF323 antibody (orb164754)

Description	Rabbit polyclonal antibody to ZNF323
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ZSCAN31
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47293
Uniprot ID	Q96LW9
NCBI	NP_001230172.1 , NP_112161.3 , NP_001128688.1 , NP_001230170.1 , NP_001230173.1 , NP_001230171.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.