

## Product Datasheet

# Recombinant Human 5'-AMP-activated protein kinase subunit beta-2(PRKAB2) (orb1647690)

|                        |                                                                                                                                                                                                                                                                                 |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Recombinant Human 5'-AMP-activated protein kinase subunit beta-2(PRKAB2)                                                                                                                                                                                                        |
| <b>Species/Host</b>    | E. coli                                                                                                                                                                                                                                                                         |
| <b>Storage</b>         | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| <b>Tag</b>             | N-terminal GST-tagged                                                                                                                                                                                                                                                           |
| <b>Note</b>            | For research use only                                                                                                                                                                                                                                                           |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                                                                                                                                                                                                  |