

Product Datasheet

BTN3A3 antibody (orb164826)

Description	Rabbit polyclonal antibody to BTN3A3 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	BTN3A3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	65002
Uniprot ID	O00478
NCBI	NP_932078.2 , NP_008925.1 , NP_001229732.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.