

Product Datasheet

RNF14 antibody (orb164838)

Description	Rabbit polyclonal antibody to RNF14 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RNF14
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	53837
Uniprot ID	Q9UBS8
NCBI	NP_899646.1 , NP_004281.1 , NP_001188294.1 , NP_899647.1 , NP_899648.1 , NP_899645.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.