

Product Datasheet

GTF2A1 antibody (orb164846)

Description	Rabbit polyclonal antibody to GTF2A1 (C-term).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GTF2A1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	41514
Uniprot ID	P52655
NCBI	NP_963889.1 , NP_056943.1 , NP_001265869.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.