


Product Datasheet

LARP1B antibody (orb164874)

Description	Rabbit polyclonal antibody to LARP1B (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	LARP1B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	105322
Uniprot ID	Q659C4
NCBI	NP_835144.1 , NP_115615.2 , NP_001265533.1 , NP_060548.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



western blot analysis of ZR-75-1 cell li...