

Product Datasheet

LONRF3 antibody (orb164876)

Description	Rabbit polyclonal antibody to LONRF3 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	LONRF3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	84490
Uniprot ID	Q496Y0
NCBI	NP_001027026.1 , NP_079054.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.