

Product Datasheet

ZNF613 antibody (orb164903)

Description	Rabbit polyclonal antibody to ZNF613 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ZNF613
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	70143
Uniprot ID	Q6PF04
NCBI	NP_079116.2 , NP_001026891.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.