

Product Datasheet

LZIC antibody (orb164913)

Description	Rabbit polyclonal antibody to LZIC
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	LZIC
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	21495
Uniprot ID	Q8WZA0
NCBI	NP_115744.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.