

Product Datasheet

MSL3 antibody (orb164914)

Description	Rabbit polyclonal antibody to MSL3 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MSL3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59824
Uniprot ID	Q8N5Y2
NCBI	NP_523353.2 , NP_523352.1 , NP_001180199.1 , NP_001269103.1 , NP_006791.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.