

Product Datasheet

NTAN1 antibody (orb164917)

Description	Rabbit polyclonal antibody to NTAN1 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	NTAN1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	34677
Uniprot ID	Q96AB6
NCBI	NP_001257696.1 , NP_775745.1 , NP_001257695.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.