

Product Datasheet

PPIL6 antibody (orb164920)

Description	Rabbit polyclonal antibody to PPIL6 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PPIL6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35228
Uniprot ID	Q8IXY8
NCBI	NP_001104768.2 , NP_775943.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.