

Product Datasheet

SR140 antibody (orb164930)

Description	Rabbit polyclonal antibody to SR140 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	U2SURP
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	118292
Uniprot ID	O15042
NCBI	NP_001073884.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.