

Product Datasheet

MBOAT2 antibody (orb164936)

Description	Rabbit polyclonal antibody to MBOAT2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MBOAT2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59527
Uniprot ID	Q6ZWT7
NCBI	NP_620154.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.