

Product Datasheet

SRGAP1 antibody (orb164941)

Description	Rabbit polyclonal antibody to SRGAP1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SRGAP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	124264
Uniprot ID	Q7Z6B7
NCBI	NP_065813.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.