

## Product Datasheet

### SPOCK3 antibody (orb164970)

|                            |                                                                                                                                                                                                                                                                             |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to SPOCK3 (N-term).                                                                                                                                                                                                                              |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                                                                      |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                                                |
| <b>Tested Applications</b> | WB                                                                                                                                                                                                                                                                          |
| <b>Immunogen</b>           | Synthetic Peptide                                                                                                                                                                                                                                                           |
| <b>Target</b>              | SPOCK3                                                                                                                                                                                                                                                                      |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.                                                                                                |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                       |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                                                  |
| <b>MW</b>                  | 49429                                                                                                                                                                                                                                                                       |
| <b>Uniprot ID</b>          | <a href="#">Q9BQ16</a>                                                                                                                                                                                                                                                      |
| <b>NCBI</b>                | <a href="#">NP_001191287.1</a> , <a href="#">NP_001238896.1</a> ,<br><a href="#">NP_001191281.1</a> , <a href="#">NP_058646.2</a> ,<br><a href="#">NP_001191284.1</a> , <a href="#">NP_001035249.1</a> ,<br><a href="#">NP_001191283.1</a> , <a href="#">NP_001191282.1</a> |
| <b>Dilution Range</b>      | WB: 1:1000                                                                                                                                                                                                                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                                             |