

Product Datasheet

SRCRB4D antibody (orb164987)

Description	Rabbit polyclonal antibody to SRCRB4D
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SSC4D
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60819
Uniprot ID	Q8WTU2
NCBI	NP_542782.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.