

Product Datasheet

ALG13 antibody (orb165016)

Description	Rabbit polyclonal antibody to ALG13
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ALG13
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	126056
Uniprot ID	Q9NP73
NCBI	NP_001093392.1, NP_001244166.1, NP_001161857.1, NP_001244163.1, NP_001244170.1, NP_001034299.3, NP_001244160.1, NP_060936.1, NP_001244159.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.