

Product Datasheet

GALNT9 antibody (orb165085)

Description	Rabbit polyclonal antibody to GALNT9
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	GALNT9
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	68359
Uniprot ID	Q9HCQ5
NCBI	NP_068580.2 , NP_001116108.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.