

Product Datasheet

IP6K3 antibody (orb165092)

Description	Rabbit polyclonal antibody to IP6K3 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	IP6K3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46417
Uniprot ID	Q96PC2
NCBI	NP_473452.2 , NP_001136355.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.