

Product Datasheet

Recombinant *Helicobacter pylori* Cytotoxicity-associated immunodominant antigen(cagA),partial (orb1651025)

Description	Recombinant <i>Helicobacter pylori</i> Cytotoxicity-associated immunodominant antigen(cagA),partial
Species/Host	Yeast
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.