

Product Datasheet

AP3M1 antibody (orb165105)

Description	Rabbit polyclonal antibody to AP3M1 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	AP3M1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46939
Uniprot ID	Q9Y2T2
NCBI	NP_036227.1 , NP_996895.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

