



Product Datasheet

AGPAT9 antibody (orb165119)

Description	Rabbit polyclonal antibody to AGPAT9
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GPAT3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	48705
Uniprot ID	Q53EU6
NCBI	NP_001243351.1 , NP_116106.2 , NP_001243350.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse heart tis...