

Product Datasheet

GLIPR1L1 antibody (orb165128)

Description	Rabbit polyclonal antibody to GLIPR1L1 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GLIPR1L1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	27151
Uniprot ID	Q6UWM5
NCBI	NP_689992.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

