

Product Datasheet

DNAJC30 antibody (orb165173)

Description	Rabbit polyclonal antibody to DNAJC30 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	DNAJC30
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	25961
Uniprot ID	Q96LL9
NCBI	NP_115693.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

