

Product Datasheet

MS4A8B antibody (orb165185)

Description	Rabbit polyclonal antibody to MS4A8B (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MS4A8
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	26290
Uniprot ID	Q9BY19
NCBI	NP_113645.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

