

Product Datasheet

N6AMT1 antibody (orb165187)

Description	Rabbit polyclonal antibody to N6AMT1 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	N6AMT1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	22957
Uniprot ID	Q9Y5N5
NCBI	NP_037372.3 , NP_877426.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

