

Product Datasheet

ITGB2 antibody (orb165193)

Description	Rabbit polyclonal antibody to ITGB2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	ITGB2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	84782
Uniprot ID	P05107
NCBI	NP_000202.2, NP_001120963.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

