

Product Datasheet

Recombinant Human Cell division cycle protein 23 homolog(CDC23),partial (orb1651989)

Description	Recombinant Human Cell division cycle protein 23 homolog(CDC23),partial
Species/Host	E. coli
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.