

Product Datasheet

Recombinant Mouse Leucine-rich repeat and immunoglobulin-like domain-containing nogo receptor-interacting protein 1(Lingo1),partial (orb1652013)

Description	Recombinant Mouse Leucine-rich repeat and immunoglobulin-like domain-containing nogo receptor-intera...
Species/Host	E. coli
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.