

## Product Datasheet

# Recombinant Human Protein disulfide-isomerase TMX3(TM3) (orb1652129)

|                        |                                                                                                                                                                                                                                                                                                           |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Recombinant Human Protein<br>disulfide-isomerase TMX3(TM3)                                                                                                                                                                                                                                                |
| <b>Species/Host</b>    | E. coli                                                                                                                                                                                                                                                                                                   |
| <b>Storage</b>         | The shelf life is related to many<br>factors, storage state, buffer<br>ingredients, storage temperature<br>and the stability of the protein<br>itself. Generally, the shelf life of<br>liquid form is 6 months at -20°C/-<br>80°C. The shelf life of lyophilized<br>form is 12 months at -20°C/-<br>80°C. |
| <b>Tag</b>             | N-terminal 6xHis-SUMO-tagged                                                                                                                                                                                                                                                                              |
| <b>Note</b>            | For research use only                                                                                                                                                                                                                                                                                     |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                                                                                                                                                                                                                            |