

Product Datasheet

ALG11 antibody (orb165215)

Description	Rabbit polyclonal antibody to ALG11 (C-term).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ALG11
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55651
Uniprot ID	Q2TAA5
NCBI	NP_001004127.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

